



PONTIFICIA
UNIVERSIDAD
CATÓLICA
DE CHILE




Instituto de Ingeniería Biológica y Médica

Scaling Up Bioimaging with Microfluidic Chips

Asst. Prof. Dr. Tobias Wenzel




Microfluidic chip - What does 'chip' mean here?




Making Chips: Photolithography and soft lithography


Photolithography



Soft lithography



iibm Instituto de Ingeniería Biológica y Médica



What Can You Do With a Chip on a Microscope?

Basic Applications



Gravity-powered systems for straightforward biological observations.

Intermediate Systems



Traps, gradients, and interface control for complex experimental conditions.

Advanced Platforms



Droplets, microgels, and organ-on-chip for cutting-edge research.



Flow chambers

Simple chambers to be mounted on advanced microscopes, and added controls

Simple Yet Powerful

- Gravity-driven flow is an option
- Ideal for observing biofilms, plant roots, and phase transitions over time
- Perfect for live cell imaging applications

VOL. 60 1994

MULTICELI






FIG. 1. Schematic diagram of the flow cell construction used to culture biofilms for microscopic examinations. The use of flow cells containing up to 10 growth chambers simplified handling and experimental replication when comparisons between biofilms grown on different carbon sources were made.




Täuber, Sarah, et al. "How to perform a microfluidic cultivation experiment—A guideline to success." *Biosensors* 11.12 (2021): 485.

Microfluidics cultivation chambers - trouble shooting






Täuber, Sarah, et al. "How to perform a microfluidic cultivation experiment—A guideline to success." *Biosensors* 11.12 (2021): 485.

Cell adhesion studies with flow control



Wei, Guanju, and Judy Q. Yang. "Microfluidic investigation of the impacts of flow fluctuations on the development of *Pseudomonas putida* biofilms." *npj Biofilms and Microbiomes* 9.1 (2023): 73.

The mother machine




Niederholtmeyer, H. et al. Rapid cell-free forward engineering of novel genetic ring oscillators. *Elife* 4, e09771 (2015).


Wang, Ping, et al. "Robust growth of *Escherichia coli*." *Current biology* 20.12 (2010): 1099-1103.

Long-term single-cell time-lapse microscopy - by type

A *B. subtilis* device with long channels and feeding side-channels




B *S. cerevisiae* device to accommodate budding in either orientation





Allard P, Papazotos F and Potvin-Trottier L (2022), Microfluidics for long-term single-cell time-lapse microscopy: Advances and applications. Front. Bioeng. Biotechnol. 10:968342. doi: 10.3389/fbioe.2022.968342

Gradients: Concentration testing in growth chambers




Yin, Bao-Sheng, et al. "An integrated microfluidic device for screening the effective concentration of locally applied tacrolimus for peripheral nerve regeneration." *Experimental and therapeutic medicine* 9.1 (2015): 154-158.



Englert DL,
Manson MD,
Jayaraman A.
Flow-based
microfluidic
device for
quantifying
bacterial
chemotaxis in
stable, competing
gradients. *Appl
Environ Microbiol.*
2009
Jul;75(13):4557-6
4. doi: 10.1128/
AEM.02952-08.

Diffusion between laminar flows - organs on a chip




Right: P. Yager et al. (2006). *Nature* 442:
412-18.

Left: P.J.A. Kenis et al. (1999). *Science*
285: 83-5.

Gas control


A Oxygen control (0% - 100% O₂) with microfluidic cultivation



B Oxygen sensing (0% - 100% O₂)


RTDP FLIM
O₂ imaging

C Automated image acquisition & analysis




Aerobic ($\mu_{aerobic} = 0.59 \text{ h}^{-1}$)

A Fluorescence images




B Generation




Anaerobic ($\mu_{anaerobic} = 0.50 \text{ h}^{-1}$)


D Fluorescence images




E Generation



C




F




Kasahara, Keitaro, et al.
"Enabling oxygen-controlled microfluidic cultures for spatiotemporal microbial single-cell analysis." *Frontiers in microbiology* 14 (2023): 1198170.

Microfluidics for precise manipulation of model organisms

a Behavioral Analysis




b Mechanical Stimulation




c Genetic Perturbation



d Chemical Stimulation



e Thermal Perturbation



Frey, N., Sönmez, U.M., Minden, J. *et al.* Microfluidics for understanding model organisms. *Nat Commun* **13**, 3195 (2022).
<https://doi.org/10.1038/s41467-022-30814-6>

- Studies of microbiome-host interactions without mice
- Greater control over sample space, time, and flexibility






Image credit: Hyun Jung Kim, Dongeun Huh, Geraldine Hamilton and Donald E. Ingber, human gut-on-a-chip inhabited by microbial flora that experiences intestinal peristalsis-like motions and flow, 2012


Organoids-on-a-chip: Modeling organ–organ interactions



Jin, Yoonhee, et al.
"Vascularized liver organoids generated using induced hepatic tissue and dynamic liver-specific microenvironment as a drug testing platform."
Advanced Functional Materials 28.37 (2018): 1801954



Trap chips

Cultivation control in chips with traps and flow








Jin, Byung-Ju, et al. "Microfluidics platform for measurement of volume changes in immobilized intestinal enteroids." *Biomicrofluidics* 8.2 (2014)

A second layer for bubble trapping (zebrafish embryos)



Zhu, Z., et al. "A Bubble-Free Microfluidic Device for Easy-to-Operate Immobilization, Culturing and Monitoring of Zebrafish Embryos. Micromachines 10 (3), 168 (2019)."

Zebrafish trapping for antibody staining cycles



Ye, Songtao, Wei-Chun Chin, and Chih-Wen Ni. "A multi-depth spiral milli fluidic device for whole mount zebrafish antibody staining." *Biomedical Microdevices* 25.3 (2023): 30.






Fig. 4. Bacterial competition for root surface colonization. Real-time im-


Droplet-Based High-Throughput Assays



Microfluidic droplets - record-keeping high-throughput




Running a microfluidic experiment (DropSeq)



Droplet Microfluidics

- Ultra-high throughput
- Single cell control
- Less contamination
- Versatile
- Reagent efficient




Droplet-based high-throughput cultivation for accurate screening of antibiotic resistant gut microbes

William J Watterson^{1,2*}, Melikhan Tanyeri^{1,2,3}, Andrea R Watson⁴,
Candace M Cham⁴, Yue Shan⁴, Eugene B Chang⁴, A Murat Eren^{4,5,6*}, Savaş Tay^{1,2*}





2020





Droplet digital PCR (ddPCR)

Applications: absolute abundance, viral load, genome copy number, etc.




- Li, H. *et al.* Application of droplet digital PCR to detect the pathogens of infectious diseases. *Bioscience Rep* **38**, BSR20181170 (2018).
- Hindson, B. J. *et al.* High-Throughput Droplet Digital PCR System for Absolute Quantitation of DNA Copy Number. *Anal Chem* **83**, 8604–8610 (2011).
- Martinez-Hernandez, F. *et al.* Droplet Digital PCR for Estimating Absolute Abundances of Widespread Pelagibacter Viruses. *Front Microbiol* **10**, 1226 (2019).
- Hindson, C. M. *et al.* Absolute quantification by droplet digital PCR versus analog real-time PCR. *Nat Methods* **10**, 1003–1005 (2013).
- Yang, R., Paparini, A., Monis, P. & Ryan, U. Comparison of next-generation droplet digital PCR (ddPCR) with quantitative PCR (qPCR) for enumeration of *Cryptosporidium* oocysts in faecal samples. *Int J Parasitol* **44**, 1105–1113 (2014)




With the highest degree of interaction, syntrophic growth can be substantially hindered in larger droplets to the extent where a significant subpopulation fails to establish co-growth .

Interactions are best studied in droplets!






Variation of: Leveraging interactions in microfluidic droplets for enhanced biotechnology screens
C. Vitalis & T. Wenzel; Current Opinion in Biotechnology, 2023/08

Many interactions can be studied in droplets!






Leveraging interactions in
microfluidic droplets for
enhanced biotechnology screens
C. Vitalis & T. Wenzel
Current Opinion in
Biotechnology, 2023/08


- High throughput analysis of host-pathogen interactions in droplets



Gel-Microdroplets and Semi-Permeable Capsules





3d-printed cellular interactions



THE BUBBLE MACHINES


Researchers use microfluidic chips to make liposomes which are similar to the envelopes that contain cellular components. The chip features a six-way junction that can fill liposomes with an inner solution. With the fatty alcohol 1-octanol in the outer solution, liposomes form around the inner solution. Over time, excess 1-octanol is removed from one end and spontaneously split off, leaving



©nature

THE ASSEMBLY LINES

A pico-injection system allows researchers to load cell-membrane-like compartments called liposomes with functional proteins. Liposomes are stabilized by a polymer coating and pushed through a microfluidic channel. As they pass over a pico-injection site, an electrical pulse can trigger the incorporation of internal proteins or membrane-bound proteins (as shown) into the liposomes.



©nature

Open source tools from our lab

Making microfluidics accessible across Latin America

Open Source Syringe Pumps and Controller

- Low-cost
- Precise low-flows (with gearbox)
- Allows manual intervention
- And remote automation



Design files, documentation and acknowledgements: https://wenzel-lab.github.io/syringe-pumps-and-controller/2_syringe-pump.html

Low-Cost Open Hardware Droplet Workstation



Strobe-enhanced microscopy stage

by Pierre Padilla-Huamantínco, Matías Hurtado-Labarca, and Tobias Wenzel
Latin American Hub for Bioimaging Through Open Hardware (LIBRE hub)

Strobe-enhanced microscopy stage

Build the 3-level microscopy stage

Print the plastic parts

Laser cut the acrylic
parts


Assemble the focus
mechanism

Assemble the basics
optics module

Attach parts to the
top plate



Attach parts to the
middle plate

3-level microscopy stage



Before you start building the station, you will need to source all the components listed in our [bill of materials](#) ([HTML](#), [CSV](#)), which is given on the next page.




Instructions



Reference: Padilla-Huamantínco, P., Durán, E., Wenzel, T. Plasmid Stability Analysis with Open-Source Droplet Microfluidics. J. Vis. Exp. (214), e67659, doi:10.3791/67659 (2024).

Open source cytometry & droplet sorting (FADS)

- Open microscope alignment stage(s)
- High-speed imaging with strobe
- Realtime FPGA electronics
- High voltage amp and electrodes on-chip
- Optomechanics
- Fluidics





Wenzel Lab

Wenzel Lab at IIBM, Universidad Católica de Chile, check also: <https://github.com/LIBREhub>

Unfollow

23 followers

Santiago, Chile

<https://wenzel-lab.github.io/en/>

@WenzelLab

wenzellab

company/92802424

<https://www.printables.com/@Wenz...>

Popular repositories

droplet-sorter-master

Public

Master repository that documents how the different parts come together of the open source hardware Fluorescence Activated (Microfluidic) Droplet Sorter (FADS).

☆ 14

modular-microfluidics-workstation-controller

Public

Repository for the development of the modular parts of a free and open source microfluidic workstation.

G-code

☆ 9

🔗 2

droplet-sorting-FPGA-controller

Public

IT setup for realtime analysis of microfluidic droplets with FPGA board (RedPitaya) based on python pyprl library

Python

☆ 8

syringe-pumps-and-controller

Public

Design of 3d printable syringe pumps and dual controller.

JavaScript

☆ 5

🔗 2

high-voltage-pulse-generator

Public

Custom electronics solution for microfluidic droplet sorting pulses

flow-microscopy-platform

Public

Our open hardware flow microscopy platform is a more comprehensive version of our strobe-enhanced microscopy stage design, including the

View as: Public

You are viewing the README and pinned repositories as a public user.

You can [create a README file](#) or [pin repositories](#) visible to anyone.

[Get started with tasks](#) that most successful organizations complete.

Discussions

Set up discussions to engage with your community!

[Turn on discussions](#)

People



We share many 3D printable designs

English

Search

244

Printables

3D Models

Store

Clubs

Community

Contests

Brands

Events

Groups

Education

+ Create

About

3D Models40

Makes25

Activity

Collections9

Lvl 18

LIBRE hub

@WenzelLab

4,317

802

87 followers

6 following

4

5

4

Advanced

Expert

White Dwarf

Follow

Message

X

Instagram

LinkedIn


GitHub

General Terms and Conditions of Use of the PRUSA Websites

Terms of Service of Printables.com

Highlighted models

All models




Adaptable Pipette Holders

62

5

261




Strobe-enhanced microscopy stage

18

5

41




Moldular tube holders

35

5

147




96 well-plate locator stand

62

4.5

351




Open-Source Syringe Pumps

71

5

142





Optical filter cubes - open-UC2 (improved)

28



0

95



Building a LatAm science maker community




Search




+ Create







Home



Shorts



Subscriptions




You

LIBRE hub

@librehub · 538 subscribers · 75 videos

The Latin American Hub for Bioimaging Through Open Hardware (LIBRE_hub) is a training ...more

librehub.github.io and 4 more links

 Subscribed


Home

Videos

Live


Playlists

Posts




Created playlists


Sort by




Taller Obtaining your first spin-echo with a low-cost ...
View full playlist




Taller MRI de bajo costo
View full playlist



Tutoriales
View full playlist



Napari for Bioimage Analysis Workshop
View full playlist



LIBRE_hub Seminars
View full playlist

Chan Zuckerberg Initiative



Instituto de Ingeniería Biológica y Médica

- Tobias Wenzel
#OpenHardware
- Pablo Irarrázaval
#openMRI
- Vicente Parot
#compMicroscopy

Latin American Hub for Bioimaging through Open Hardware (LIBRE hub)

